

ERYTHROCYTE
TOTAL PANEL
(Total,
Essential and
Toxic)

Orderable – TEESSERC (Essential Panel) TETOXERC (Toxic Panel) TEERC (Essential and Toxic panel)

Turn Around Time: 10 days

Specimen:

6 mL K2-EDTA Royal Blue Vacutainer tube

Collection Information:

[Trace Elements Collection Requirements](#)

Reference Ranges:

[Find individual Reference Ranges here:](#)

Reference Ranges are based on Non-Occupationally exposed population.

Interpretive Comments:

Find Interpretive Comment and Clinical Information [here:](#)



Laboratory:
Trace Elements Lab



Requisition:
TRACE ELEMENTS
REQUISITION



Method of Analysis:
High Resolution Sector
Field Inductively Coupled
Mass Spectrometry
[HR-SF-ICP-MS](#)



Test Schedule:
Batched Analysis