

# BERYLLIUM, URINE

## Orderable - BEU / BEU24

Turnaround Time: 10 days

### Specimen:

24-hour urine collected in an unused 24-hour urine container or random urine

### Collection Information:

#### [Trace Elements Collection Requirements](#)

### Reference Ranges:

#### SI Units:

Random Urine:  
0.0-2.4 umol/mol creatinine

24 Hour Urine:  
0-16.7 nmol/d

#### Conventional Units:

Random Urine:  
0-191 ng/g creatinine

24 Hour Urine:  
0-0.150 µg/d

Reference Ranges are based on Non-Occupationally exposed population.

### Interpretive Comments:

Find Interpretive Comment and Clinical Information [here:](#)

### Specimen Stability:

Ambient: 14 days  
Refrigerated: 11 months  
Frozen: 11 months



**Laboratory:**  
Trace Elements Lab



**Requisition:**  
TRACE ELEMENTS  
REQUISITION



**Method of Analysis:**  
ICP-MS/MS



**Test Schedule:**  
Batched analysis