

CALCIUM, WHOLE BLOOD

(tested in Trace
Elements)

Orderable - CABLD

Turnaround Time: 10 days

Specimen:

6 mL K2-EDTA Royal Blue Vacutainer tube

Collection Information:

[Trace Elements Collection Requirements](#)

Reference Ranges:

SI Units (Reported on Patient Chart):	Conventional Units:
1.06-1.29 mmol/L	42.6-51.8 mg/L

Reference Ranges are based on Non-Occupationally exposed population.

Interpretive Comments:

Find Interpretive Comment and Clinical Information [here](#):

Specimen Stability:

Ambient: 22 days
Refrigerated: 14 months
Frozen: 14 months



Laboratory:
Trace Elements Lab



Requisition:
TRACE ELEMENTS
REQUISITION



Method of Analysis:
ICP-MS/MS



Test Schedule:
Batched analysis