

MERCURY, URINE

Orderable – HGU / HGU24

Turnaround Time: 10 days

Specimen:

24-hour urine collected in new 24-hour urine container or random urine

Collection Information:

[Trace Elements Collection Requirements](#)

Reference Ranges:

SI Units:

Random Urine:
0-1.2 $\mu\text{mol}/\text{mol}$ creatinine

24 Hour Urine:
0-20.0 nmol/d

Conventional Units:

Random Urine:
0-2.1 $\mu\text{g}/\text{g}$ creatinine

24 Hour Urine:
0-4.0 $\mu\text{g}/\text{d}$

Reference Ranges are based on Non-Occupationally exposed population.

Interpretive Comments:

Find Interpretive Comment and Clinical Information [here](#):

Specimen Stability:

Ambient: 14 days
Refrigerated: 11 months
Frozen: 11 months



Laboratory:
Trace Elements Lab



Requisition:
TRACE ELEMENTS
REQUISITION



Method of Analysis:
ICP-MS/MS



Test Schedule:
Batched analysis