

THALLIUM,  
ERYTHROCYTES

## Orderable - TLERC

Turn Around Time: 10 days

### Specimen:

6 mL K2 EDTA Royal Blue Vacutainer tube

### Collection Information:

[Trace Elements Collection Requirements](#)

### Reference Ranges:

| SI Units (Reported on Patient Chart): | Conventional Units: |
|---------------------------------------|---------------------|
| 0.00-0.13 nmol/L                      | 0.000-0.026 µg/L    |

Reference Ranges are based on Non-Occupationally exposed population.

### Interpretive Comments:

Find Interpretive Comment and Clinical Information [here](#):

### Specimen Stability:

Ambient: 22 days  
Refrigerated: 14 months  
Frozen: 14 months



**Laboratory:**  
Trace Elements Lab



**Requisition:**  
TRACE ELEMENTS  
REQUISITION



**Method of Analysis:**  
High Resolution Sector  
Field Inductively Coupled  
Mass Spectrometry  
[HR-SF-ICP-MS](#)



**Test Schedule:**  
Batched analysis