

VANADIUM,  
ERYTHROCYTES

## Orderable - VERC

Turn Around Time: 10 days

### Specimen:

6 mL K2 EDTA Royal Blue Vacutainer tube

### Collection Information:

[Trace Elements Collection Requirements](#)

### Reference Ranges:

SI Units (Reported on Patient Chart):	Conventional Units:
0.2-1.1 nmol/L	0.012-0.054 µg/L

Reference Ranges are based on Non-Occupationally exposed population.

### Interpretive Comments:

Find Interpretive Comment and Clinical Information [here](#):

### Specimen Stability:

Ambient: 22 days  
Refrigerated: 14 months  
Frozen: 14 months



**Laboratory:**  
Trace Elements Lab



**Requisition:**  
TRACE ELEMENTS  
REQUISITION



**Method of Analysis:**  
High Resolution Sector  
Field Inductively Coupled  
Mass Spectrometry  
[HR-SF-ICP-MS](#)



**Test Schedule:**  
Batched analysis