

ZINC  
PROTOPORPHYRIN,  
WHOLE  
BLOOD

## Orderable - ZPP

Turn Around Time: 10 days

### Specimen:

6 mL K2 EDTA Royal Blue Vacutainer tube

### Collection Information:

[Trace Elements Collection Requirements](#)

### Reference Ranges:

SI Units (Reported on Patient Chart):	Conventional Units:
0.0-20.9 nmol/L	0.0-1.0 µg/L

Reference Ranges are based on Non-Occupationally exposed population.

### Interpretive Comments:

Find Interpretive Comment and Clinical Information [here](#):



**Laboratory:**  
Trace Elements Lab



**Requisition:**  
TRACE ELEMENTS  
REQUISITION



**Method of Analysis:**  
Proto Fluor



**Test Schedule:**  
Batched analysis