

AMOEBIC SEROLOGY

Orderable – PARASRO

Turnaround Time: Up to 10 days

Specimen Type and Collection Information:

Please review the Public Health Ontario Laboratory Test Information Index:

Amoebic - Serology

Specimen Type	Minimum Volume	Collection Kit
Serum	1.0 mL	Vacutainer tubes (SST) - Gold 5mL
Blood, clotted	5.0 mL	Vacutainer tubes (SST)- Gold 5mL

An acute (collected early after the onset of symptoms) and a convalescent (collected 2-3 weeks later) may be required for laboratory diagnosis.

Critical Information Required:

Specify suspected parasite in the **Specimen Description** field when placing order in PowerChart and on the Public Health Laboratory Requisition.

Storage and Shipment:

The specimen must be refrigerated.
The specimen must be received at Microbiology within 24 hours of collection.



Laboratory:
Microbiology



Requisition:
[PUBLIC HEALTH
LABORATORY TEST
REQUISITION](#)

This requisition **MUST**
accompany the specimen



Method of Analysis:
EIA



Test Schedule:
Referred out weekdays



Referred to:
Public Health Laboratory
1200 Commissioners Rd E
London, Ontario
519-455-9310