

HTLV I and II SEROLOGY

Orderable – HTLV₁, HTLV₂

Turnaround Time: Up to 14 days

Alternate Name(s):

Human T-cell Lymphotropic Virus

Specimen Type and Collection Information:

Please review the Public Health Ontario Laboratory Test Information Index:

[HTLV I and II - Serology](#)

Specimen Type	Minimum Volume	Collection Kit
Serum	1.0 mL	Gold/Red Vacutainer tubes (SST)
Blood	3.0 mL	Gold/Red – vacutainer tubes (SST)

Storage and Shipment:

The specimen must be refrigerated.
The specimen must be received within 24 hours of collection.



Laboratory:
Microbiology



Requisition:
[PUBLIC HEALTH
LABORATORY TEST
REQUISITION](#)

This requisition **MUST**
accompany the specimen



Method of Analysis:
CMIA
LIA



Test Schedule:
Referred out weekdays



Referred to:
Public Health Laboratory
1200 Commissioners Rd E
London, Ontario
519-455-9310