

HISTOPLASMA ANTIGEN

Orderable – HISTOAG

Turnaround Time: Up to 6 days

Specimen and Collection Information:

Specimen Type	Minimum Volume	Collection Kit
Urine	0.5 mL	Sterile Container



Laboratory:
Microbiology



Requisition:
[MICROBIOLOGY
Requisition](#)



Method of Analysis:
Quantitative Sandwich
Enzyme Immunoassay
(EIA)



Test Schedule:
Referred out weekdays



Referred to:
MiraVista Diagnostics
4705 Decatur Blvd.
Indianapolis, Indiana
46241
USA

Interpretive Comments:

Positive: 0.20-20:00 ng/mL

Positive Above the Limit of Quantification:

Results greater than 20.0 ng/mL fall outside the linear range of the assay. These results are positive, but not accurately quantifiable.

Limitations

- Result is not intended to be used as the sole means for clinical diagnosis or patient management decisions
- Cross-reactions are seen with blastomycosis, paracoccidioidomycosis, penicilliosis, less frequently in coccidioidomycosis, rarely in aspergillosis and possibly sporotrichosis

Comments:

Testing restricted to Infectious Diseases (ID) service following a Microbiologist consultation.

Storage and Shipment:

Refrigerate (store), Room temperature (ship)