

HLA TRANSFUSION TYPING

Orderable - HLA Transfusion Typing (ABC)

Turnaround Time: 1-3 days

STAT: Same day if blood is received in the laboratory before 0900

Alternate Name(s):

Transfusion ABC

Specimen:

Adult	Pediatric
2 x 4 mL Lavender (EDTA) top Vacutainer tubes	<1 year: 2 mL EDTA 1-10 years: 2 mL EDTA 10-18 years: 4 mL EDTA
Note: Increase collection volumes if leukocyte count is low, smaller volumes apply for paediatric patients	Contact the laboratory if these volumes are not possible: 519-663-3320

Collection Information:

Samples that are being shipped into the laboratory can be drawn Monday- Friday

Draw samples prior to treatment if chemotherapy is being considered.

Reference Ranges:

See report

Comments:

This test is performed in preparation for receiving HLA matched platelets.

Results will be faxed upon request.



Laboratory:
Transplant Lab



Requisition:
[TRANSPLANT
LABORATORY
REQUISITION](#)



Method of Analysis:
HLA Class I typings are performed using molecular techniques (RSSO or SSP).



Test Schedule:
As required Monday to Friday except Stat holidays

A large graphic on the left side of the page consists of a dark blue outer ring and a light blue inner circle. Inside the inner circle, the text "HLA TRANSFUSION TYPING" is written in white, bold, sans-serif capital letters.

HLA TRANSFUSION TYPING

Critical Information Required:

Indicate the purpose of the test.

Storage and Shipment:

Keep EDTA blood at room temperature and deliver within 72 hours.