

FARMER'S
LUNG IGG
ANTIBODIES
SERUM

Test currently unavailable due to reagent shortage (2022/05/05)

Orderable - AAFL

Turn Around Time: 10 days

Alternate Name(s):

Allergic Alveolitis
Allergic Lung Serology
Farmer's Lung
Farmers Lung
Farmer's Lung Precipitins
Micropolyspora faeni
M. faeni
Thermoactinomyces vulgaris
T. vulgaris
Hypersensitivity Pneumonitis



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Fluorescence Enzyme
Immunoassay (FEIA)



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	2-10 years: 3 mL Red top Vacutainer tube


Collection Information:

A single 1 mL aliquot is sufficient to perform this test as well as Aspergillus fumigatus, Avian, Budgie, Pigeon and Parrot IgG Antibodies tests.

To order Allergic Lung Serology please order both Farmer's Lung IgG Antibodies AND Aspergillus fumigatus IgG Antibodies (separate orderable).

Reference Ranges:

M. faeni IgG antibodies:	<5.2 mg/L
T. vulgaris IgG antibodies:	<21.5 mg/L



FARMER'S
LUNG IGG
ANTIBODIES
SERUM

Storage and Shipment:

Store and ship and serum frozen.

Specimen Stability:

Refrigerated: 7 days

Frozen: 3 months