

HLA B\*15:02 CBZ  
HYPERSENSITIVITY  
CASE

# Orderable - HLA B\*15:02 CBZ Hypersensitivity case

Turnaround Time: 14 days

## Specimen:

2 x 4 mL Lavender (EDTA) top Vacutainer tubes

## Collection Information:

Samples that are being shipped into the laboratory can be drawn Monday- Friday

## Reference Ranges:

See report

## Storage and Shipment:

Keep EDTA blood at room temperature and deliver within 72 hours.



**Laboratory:**  
Transplant Lab



**Requisition:**  
[TRANSPLANT  
LABORATORY  
REQUISITION](#)



**Method of Analysis:**  
HLA Class I and Class II  
 typings are performed  
 using molecular  
 techniques (RSSO or SSP)



**Test Schedule:**  
As required Monday to  
Friday except Stat  
holidays