

ADENOVIRUS QUANTITATIVE PCR

Orderable - ADENOQ

Turnaround Time: 4 business days from receipt at testing lab

Specimen Type and Collection Information:

Specimen Type	Minimum Volume	Collection Kit
Plasma	500 uL	EDTA (purple top vacutainer)



Laboratory:
Microbiology



Requisition:
Electronic order



Method of Analysis:
PCR



Test Schedule:
Referred out weekdays



Referred to:
Sick Kids
555 University Avenue
Room3676, Atrium
Toronto, ON M5G 1X8
416-813-7200

Special Processing:

Adenovirus PCR test is available for bone marrow transplant patients only. Any other request requires approval by the Microbiologist on-call (Pager #19000).

Storage and Shipment:

Blood (EDTA): Tubes MUST be spun, separated and stored at 4°C for shipment the next morning if there is a delay of ≥4 hours in transport to the Microbiology Laboratory. If there is a delay of more than 24 hours, plasma MUST be frozen at -80°C.