



MITOGEN STIMULATION ASSAY

Orderable - LYMPRO

Alternate Name(s):

Lymphocyte Mitogen
Lymphocyte Proliferation Assay

Specimen:

2 X 6 mL Dark Green top Vacutainer tubes (Sodium Heparin)

Collection Information:

Sample must be drawn Wednesday.

Unauthorized samples will be discarded.

Reference Ranges:

Included in report

Special Processing:

Sample must arrive by 9 am next day delivery.

Shipping details (tracking number) must be sent to Dr. Ben Hedley for tracking purposes.



Laboratory:
Core Lab



Requisition:
Available upon request.
Request and
authorization to Dr. Ben
Hedley (ext. 55402)



Test Schedule:
Weekly (authorized
samples must be in Lab
by 10 am Wednesday).
Unauthorized samples
will be discarded.



Referred Out Location:
Alberta Precision
Laboratories