

ANTIENDOMYSIAL ANTIBODIES

Orderable - EMA

Turnaround Time: 25 days

Alternate Name(s):

EMA
Anti Mycelial

Specimen:

6 mL Red top Vacutainer tube

Collection Information:

A minimum volume of 1 mL of serum is required for testing.

Note:

This test is not to be used for celiac disease screening unless the patient is TTGAb negative and IgA deficient. **This test is approved for Dr. Andreanne Zizzo and Dr. Dhandapani Ashok.** All others will require biochemist approval to order this test; approval can be requested by calling the clinical biochemist at extension 77626 or the biochemist on-call.

Reference Ranges:

Negative

Specimen Stability:

Refrigerated: 7 days
Frozen: >7 days



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Immunofluorescence
Assay (IFA)



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
In-Common Laboratories