

DOXEPIN +
METABOLITE,
SERUM/
PLASMA

Orderable - DOX

Turnaround Time: 5 days

Alternate Name(s):

Sinequan

Specimen:

6 mL Red top Vacutainer tube

Lavender or Pink (EDTA) top tubes are also acceptable

**Avoid blood collection tubes with separator gels (e.g., Gold top).
Light Green (Li-Heparin) top tubes are also NOT acceptable.**

Collection Information:

Minimum volume required is 3 mL.

To monitor therapy, draw trough specimen prior to morning dose or 10 - 12 hours after last drug administration. Separate and freeze as soon as possible.

Reference Ranges:

Doxepin + Desmethyldoxepin:	50-150 ng/mL
--------------------------------	--------------

Critical Value:

Toxic:

Doxepin + Desmethyldoxepin:	> 300 ng/mL
--------------------------------	-------------

Storage and Shipment:

Send frozen samples ASAP.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
LC-MS/MS



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)