

BARTONELLA SEROLOGY

Orderable - CATSF

Turnaround Time: Up to 10 days

Alternate Name(s):

Cat Scratch Fever (Bartonella) Antibody

Specimen Type and Collection Information:

Please review the Public Health Ontario Laboratory Test Information Index:

[Bartonella henselae](#)

Specimen Type	Minimum Volume	Collection Kit
Serum	1.0 mL	Vacutainer tubes (SST), Gold 5 mL
Blood, Clotted	5.0 mL	Vacutainer tubes (SST), Gold 5 mL

Storage and Shipment:

The specimen must be refrigerated.
The specimen must be received within 24 hours of collection.



Laboratory:
Microbiology



Requisition:
See electronic order



Method of Analysis:
IFA IgG



Test Schedule:
Referred out weekdays



Referred to:
Public Health Laboratory
1200 Commissioners Rd E
London, Ontario
519-455-9310