

CHORIONIC
GONADOTROPIN,
FLUID

Orderable - HCGF

Turnaround Time: 4 hours

Alternate Name(s):

hCG
Beta hCG

Specimen:

Fluid

Collection Information:

Collect a minimum of 1 mL of fluid in a sterile tube without additive (e.g. orange top urine container). Grossly hemolyzed samples are unacceptable.

May also wish to collect a Light Green top (Li-Heparin) Vacutainer tube of blood to compare fluid hCG result to plasma hCG result.

Reference Ranges:

No reference range available for fluid

Interpretive Comments:

It may be helpful to interpret fluid hCG result relative to plasma hCG result.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Roche
Electrochemiluminescence



Test Schedule:
As required