

ALDOSTERONE,
URINE

Orderable - ALDU24

Turn Around Time: 2 weeks

Specimen:

24-hour urine

Collection Information:

[Collection of 24-Hour Urine Sample](#)

Collect urine over a 24-hour period in a container with no preservative.

Patient must be on a diet with normal sodium intake and not receiving diuretics for one week before urine sample is collected for this test.

Reference Ranges:

≥14 years: 3-78 nmol/day

Comments:

Aldosterone output increases in normal persons on low salt intake or on diuretic medication.

Critical Information Required:

Transfer a 50-mL aliquot of the measured urine to a sterile container. Indicate the 24 hour total volume collected.

Storage and Shipment:

Store and ship frozen.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Diasorin
Chemiluminescence



Test Schedule:
Referred out Monday-
Friday



Referred Out Location:
Dynacare, London