

ALPHA  
FETOPROTEIN  
(NON-  
PREGNANCY),  
FLUID

## Orderable - AFPF

Turn Around Time: 4 hours

### Alternate Name(s):

Alphafetoprotein  
AFP

### Specimen:

Fluid

### Collection Information:

Collect a minimum of 1 mL of fluid and transfer into a Red top Vacutainer tube. Grossly hemolyzed samples are unacceptable.

May also wish to collect a Light Green top (Li-Heparin) Vacutainer tube of blood to compare fluid AFP result to plasma AFP result.

### Reference Ranges:

No reference range available for fluid

### Interpretive Comments:

It may be helpful to interpret fluid AFP result relative to plasma AFP result.



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Roche  
Electrochemiluminescence



**Test Schedule:**  
As required