

ANTI ISLET ANTIBODIES, SERUM

Orderable - AICA

Turn Around Time: 5 days

Alternate Name(s):

Islet Cell Antibodies
Pancreatic Antibodies

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	0-2 years: 0.5 mL Gold top Microtainer 2-10 years: 2 mL Red top Vacutainer tube

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Separate serum ASAP.

Reference Ranges:

Negative

Special Processing:

Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance and should be documented.

Storage and Shipment:

Intercampus - Ship via routine or stat specimen transport. Serum can be stored for up to 4 days at 2-8°C. For longer storage, store at -20°C.

Do not ship after 1500 h.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Indirect
immunofluorescence



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[Dynacare Medical
Laboratories](#)

ANTI ISLET
ANTIBODIES,
SERUM

Store at 2-8° and ship ASAP (Monday to Friday only).

Intercity - Freeze and store serum samples at -20°C. Ship on ice packs. Avoid freeze/thaw cycles.