

ANTI-
THYROID
PEROXIDASE
ANTIBODIES,
SERUM/PLASMA

Orderable - ATPO

Turn Around Time: 5 days

Alternate Name(s):

Anti-Thyroperoxidase Antibodies
Anti-TPO Antibodies
Anti-Microsomal Antibodies

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	0-2 years: 0.5 mL Red or Gold top Microtainer 2-10 years: 3 mL Red top Vacutainer tube
Light Green top (Li-Heparin) tubes are also acceptable	
Lavender top (EDTA) tubes are not acceptable	

Collection Information:

Minimum volume of serum or plasma required is 700 µL for adult samples or 200 µL for pediatric samples.

Reference Ranges:

≤ 34 IU/mL

Comments:

Biotin may interfere with this test. Samples should not be taken from patients receiving high biotin doses (i.e. > 5 mg/day) until at least 8 hours after the last biotin administration.

Falsely elevated results may be obtained in patients being treated with Itraconazole.

Storage and Shipment:



Laboratory:
Core Lab




Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Roche
Electrochemiluminescence



Test Schedule:
Monday - Friday
0800-1600



ANTI-
THYROID
PEROXIDASE
ANTIBODIES,
SERUM/PLASMA

Store and ship serum or plasma frozen.