

ANTIENDOMYSIAL
ANTIBODIES

Orderable - EMA

Turn Around Time: 25 days

Alternate Name(s):

EMA
Anti Mycelial

Specimen:

6 mL Red top Vacutainer tube

Collection Information:

A minimum volume of 1 mL of serum is required for testing.

Note:

This test is not to be used for Celiac screening unless TTGAb negative and patient IgA deficient - Call the clinical biochemist at 77626 or biochemist on call, Monday to Friday 0800 to 1600 h for approval before ordering.

Reference Ranges:

Negative

Specimen Stability:

Refrigerated: 7 days
Frozen: >7 days**Laboratory:**
Core Lab**Requisition:**
GENERAL LABORATORY
REQUISITION**Method of Analysis:**
Immunofluorescence
Assay (IFA)**Test Schedule:**
Referred out Monday-
Thursday**Referred Out Location:**
[In-Common Laboratories](#)