

 BARBITURATE  
SCREEN

## Orderable - BARBSU

Turn Around Time: 24 hours

STAT: 1 hour

### Alternate Name(s):

Barb Screen

### Specimen:

Random urine

### Collection Information:

Collect random urine in a sterile container.

### Reference Ranges:

Negative

### Interpretive Comments:

A positive result for the Barbiturate Screen indicates that the drug may be present in urine at a minimum concentration of 200 ng/mL, and is to be regarded as preliminary/presumptive unless subsequently verified by a definitive confirmation test (available upon request by calling the Toxicology Laboratory).

**Laboratory:**  
Core Lab**Requisition:**  
GENERAL LABORATORY  
REQUISITION**Method of Analysis:**  
Kinetic Interaction of  
Particles in Solution  
(KIMS)**Test Schedule:**  
As required