

CEA, FLUID

Orderable - CEAF

Turn Around Time: 4 hours

Alternate Name(s):

Carcinoembryonic Antigen

Specimen:

Fluid

Collection Information:

Collect a minimum of 1 mL of fluid and transfer into a Red top Vacutainer tube. Grossly hemolyzed samples are unacceptable.

May also wish to collect a Light Green top (Li-Heparin) Vacutainer tube of blood to compare fluid CEA result to plasma CEA result

Reference Ranges:

No reference range available for fluid

Interpretive Comments:

It may be helpful to interpret fluid CEA result relative to plasma CEA result.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Roche
Electrochemiluminescence



Test Schedule:
As required