

COCAINE SCREEN, URINE

Orderable - COCASU

Turn Around Time: 24 hours

STAT: 1 hour

Alternate Name(s):

Cocaine

Specimen:

Random urine

Collection Information:

Collect random urine in a sterile container.

Reference Ranges:

Negative

Interpretive Comments:

A positive result for the Cocaine Screen indicates that the drug may be present in urine at a minimum concentration of 300 ng/mL, and is to be regarded as preliminary/presumptive unless subsequently verified by a definitive confirmation test (available upon request by calling the Toxicology Laboratory).



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Enzymatic immunoassay
(EIA)



Test Schedule:
As required