

# CORTISOL, 24-HOUR URINE

## Orderable – CORU24

Turn Around Time: 7 days

### Specimen:

24-hour urine

### Collection Information:

#### [Collection of 24-Hour Urine Sample](#)

Collect urine over a 24-hour period in a container with no preservative.

Please complete the laboratory label attached to the collection bottle to include patient's name, collection start time and date and collection end time and date.

### Reference Ranges:

≤ 274 nmol/day

### Special Processing:

Transfer a 10-mL aliquot of the measured urine to a sterile container and freeze.



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
IMMULITE



**Test Schedule:**  
Referred out Monday-  
Friday



**Referred Out Location:**  
[LifeLabs](#)