



CYTOMEGALOVIRUS  
SEROLOGY

## Orderable - CMVIGG/CMVIGM

Turn Around Time: 48 hours

---

### Alternate Name(s):

CMV IgM/IgG

### Specimen:

5 mL Gold or 6 mL Red top Vacutainer tube

### Collection Information:

Collect blood aseptically in a Vacutainer tube.

If congenital infection is being considered, send a sample from both the neonate and mother.

### Reference Ranges:

See report

### Interpretive Comments:

A single serum sample will yield significant results for establishing previous exposure.

### Critical Information Required:

Provide clinical history.

Indicate purpose of the request as diagnosis of acute disease (IgM) or determining previous exposure (IgG).



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Electrochemiluminescence  
Immunoassay



**Test Schedule:**  
24/7