

DOXEPINE +  
METABOLITE,  
SERUM/  
PLASMA

## Orderable - DOX

Turn Around Time: 5 days

### Alternate Name(s):

Sinequan

### Specimen:

2 x 6 mL Gold top Vacutainer or 2 x 4 mL Lavender top EDTA Vacutainer tube  
Avoid gel-separator tubes.

### Collection Information:

Minimum volume of serum required is 3 mL.

To monitor therapy, draw trough specimen prior to morning dose or 10 - 12 hours after last drug administration. Separate and freeze as soon as possible. Analysis includes Desmethyldoxepin.

### Reference Ranges:

<b>Doxepin + Desmethyldoxepin:</b>
350-1100 nmol/L

**Toxic:** Doxepin + Desmethyldoxepin: Greater than 1400 nmol/L

### Storage and Shipment:

Send frozen samples ASAP.



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
GLC/NP Detector



**Test Schedule:**  
Referred out Monday-  
Thursday



**Referred Out Location:**  
[In-Common Laboratories](#)