

ETHOSUXIMIDE
SERUM

Orderable - ETHOS

Turn Around Time: 7 days

Alternate Name(s):

Anticonvulsant
Zarontin

Specimen:

5 mL Gold or 6 mL Red top Vacutainer tube

Collection Information:

Minimum volume required is 1 mL of serum

Due to long half-life of drug time of draw is not important. When making comparative therapeutic measurements timing should be consistent.

Reference Ranges:

Therapeutic	Critical
283- 708 µmol/L	>1060 µmol/L

Critical Information Required:

Specify time of last dose and collection.

Storage and Shipment:

Store and ship refrigerated.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
High performance liquid
chromatography (HPLC)



Test Schedule:
Referred out Monday-
Friday



Referred Out Location:
[Dynacare Laboratory](#)
245 Pall Mall Street
London, Ontario
Tel: 519-679-1630