

FLUOXETINE

Orderable - FLUOX

Turn Around Time: 7 days

Alternate Name(s):

Prozac
Norfluoxetine

Specimen:

Adult	Pediatric
6 mL Red top Vacutainer tube or 4.5 mL Lavender top tube	0-2 years: 2 mL Red top
Avoid gel-separator tubes	

Collection Information:

Minimum volume required for testing is 1 mL.

Collect blood aseptically in a Vacutainer tube.

Submit trough specimen (i.e. collected within 1 hour prior to next dose.)

Assay includes the metabolite Norfluoxetine.

Reference Ranges:

Fluoxetine	Norfluoxetine
160-1600 nmol/L	170-1700 nmol/L



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
HPLC



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)