

HEPATITIS
A
ANTIBODY
IGM

Orderable - HpAIgMAb

Turn Around Time: 48 hours

Alternate Name(s):

anti HAV IgM
Hepatitis anti-HAV IgM

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer	0-2 years: 0.5 mL Gold top Microtainer 2-10 years: 2 mL Gold top tube

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Reference Ranges:

See report

Interpretive Comments:

Sample will be tested for Hepatitis A IgM antibody. Negative result indicates serology is not compatible with acute Hepatitis A infection. Positive result indicates serology is compatible with recent Hepatitis A infection.

Comments:

Heparinized samples should be avoided as they cause erroneous results.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Electrochemiluminescence
Immunoassay



Test Schedule:
24/7