



HPV TESTING,
GYNAECOLOGICAL
LIQUID BASED
PAP TEST
(CYTOLOGY)

Orderable – E-order/Requisition

Turnaround Time: 2 weeks

Alternate Name(s):

Liquid based PAP test (Cytology)
PAP test
LBP



Laboratory:
Cytopathology Lab



Requisition:
[CYTOPATHOLOGY
REQUISITION- NON-
GYNAECOLOGICAL AREA](#)



Method of Analysis:

Specimen:

Cervical, endocervical, vaginal samples

Collection Information:

Note: This test must be paid for by the patient. If the patient chooses to have this testing completed, they will be invoiced by LifeLabs directly.

Sample must be submitted in GYN PreservCyt fixative

Click here for the [External Client Cytology Collection Method](#) if you are not at LHSC or SJH

Reference Ranges:

See report

Interpretive Comments:

Clinical history is an important component for diagnostic interpretation

Special Processing:

The specimen is ThinPrep processed and must be collected in specific Gyn collection containers (PreservCyt)



Test Schedule:
Specimens are sent to Life Labs weekly for testing (Wednesday morning)

HPV TESTING,
GYNAECOLOGICAL
LIQUID BASED
PAP TEST
(CYTOLOGY)

Critical Information Required:

Specimen container must be labelled with the patient's name.

Unlabeled specimen containers cannot be processed.

Order must be accompanied by signed Cytology & HPV Testing Requisition

Storage and Shipment:

Transport all specimens in a biohazard transport bag with the specimen in the zipped area, separating the requisition in the side pocket