ERYTHROPOIETIN SERUM

Orderable - EPO

Turn Around Time: 7 days

Alternate Name(s):

EPO

Specimen:

Adult	Pediatric
5 mL Gold top	0-10 years: 2 mL Gold top tube
Vacutainer tube	

Collection Information:

Due to diurnal variation, it is recommended that specimens be collected between 0730 and noon.

Minimum volume of serum required for testing is 0.5 mL.

Plasma is not suitable.



Test Schedule: Referred out Monday-Thursday



Referred Out Location: In-Common Laboratories

Male		
Age	Range	
<1 year:	Not established	
1-3 years:	1.7-17.9 mIU/mL	
4-6 years:	3.5-21.9 mIU/mL	
7-9 years:	1.0-13.5 mIU/mL	
10-12 years:	1.0-14.0 mIU/mL	
13-15 years:	2.2-14.4 mIU/mL	
16-18 years:	1.5-15.2 mIU/mL	
>18 years:	2.6-18.5 mIU/mL	

Female		
Age	Range	



Pathology and Laboratory Medicine



Laboratory: Core Lab



Requisition: GENERAL LABORATORY REQUISITION

Method of Analysis: Immunoassay (IA)



Pathology and Laboratory Medicine

	<1 year:	Not established
	1-3 years:	2.1-15.9 mIU/mL
	4-6 years:	2.9-8.5 mIU/mL
1	7-9 years:	2.1-8.2 mIU/mL
Ĺ	10-12 years:	1.1-9.1 mIU/mL
	13-15 years:	3.8-20.5 mIU/mL
	16-18 years:	2.0-14.2 mIU/mL
	>18 years:	2.6-18.5 mIU/mL

Storage and Stability:

Stability: 10 days at room temperature, 18 days refrigerated, 28 days frozen

