

Available at VH only

London Health Sciences Centre



Pathology and Laboratory Medicine

Orderable: SSPC

Turnaround Time: 4 hours STAT: 1 hour

Alternate Name(s):

PCT

SSPC

PC



Laboratory: Core Lab



Requisition:GENERAL LABORATORY
REQUISITION



Method of Analysis: Electrochemiluminescence



Test Schedule: As required

Specimen:

Pediatric patients less than 60 days old:

0.5 mL Green top Heparin Microtainer tube

Other acceptable matrix includes serum, heparinized and EDTA plasma. Plasma/serum tubes containing separating gel can be used.

Collection Information:

Collect blood aseptically in microtainer or vacutainer tubes as per standard operating procedures.

Reference Range:

≤0.5 μg/L

Interpretive Comments:

Values ≤0.5 µg/L represent a low risk of severe bacterial infection.

Comments:

Procalcitonin testing is currently available for investigation of infants aged 0-60 days with a suspicion of a serious bacterial infection. This test **is not offered** for patients >60 days old.

Storage and Shipment:

24 hours at 20-25 °C, 48 hours at 2-8 °C, 24 months at -20 °C (\pm 5 °C). Freeze only once.