TRYPEPTIDYL PEPTIDASE 1, DRIED BLOOD SPOT/ FIBROBLAST

Orderable - DBS: TPPI/ Fibroblasts: Use Requisition to order

Alternate Name(s):

CLN2 Peptidase LINCL Batten Disease



Laboratory: Biochemical Genetics Lab

Specimen:

|||| × – Requisition: BIOCHEMICAL GENETICS LAB REQUISITION 1. Blood dot on PKU-type card

2. Fibroblast

Collection Information:

1. Saturate each circle completely (blood dot should be visible on both front and back of card) with a single application; air dry thoroughly.

2. Collect fibroblasts and send in a sterile container.



Test Schedule: As required

Reference Ranges:

Dried Blood Spot:	2.6-7.8 nmol/24 hr/3 mm disc
Fibroblast:	47-114 nmol/hr/mg protein



Pathology and Laboratory Medicine