

VALPROIC
ACID-
FREE,
SERUM

Orderable - FVAL

Turnaround Time: 5 days

Alternate Name(s):

Depakene

Specimen:

Adult	Pediatric
6 mL Red top Vacutainer tube	0-2 years: 0.5 mL Red top Microtainer 2-10 years: 4 mL Red top Vacutainer tube
Green and gold top Vacutainers with separator gels are not acceptable.	



Laboratory:
Toxicology/Special
Chemistry Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
LC-MS/MS



Test Schedule:
Monday--Friday
0800-1600

Collection Information:

Minimum volume of serum required for testing is 2 mL.

Spin and aliquot serum and freeze immediately and send to the lab.

Reference Ranges:

35.0-140.0 $\mu\text{mol/L}$

Storage and Shipment:

Separate and freeze. Ship on ice packs or frozen.