



Pathology and Laboratory Medicine

VOLTAGE
GATED
CALCIUM
CHANNEL
ANTIBODIES

Orderable – GENX: then specify test

Turnaround Time: 20 days

Alternate Name(s):

Anti VGCC

Specimen:

Adult	Pediatric
2 x 5 mL Gold top	2-10 years: 1-3 x 3.5 mL Gold top
Vacutainer tube	Vacutainer tube

Laboratory: Core Lab

Requisition:

UBC Neuro-Immunology Laboratory Requisition must be completed by neurologist



Method of Analysis: RIA

Collection Information:

2-5 mL serum required for testing

Reference Ranges:

Negative	Positive
≤30 pmol/L	>30 pmol/L

鱑

Test Schedule:Referred out Mondays

Special Processing:

Referred Out Location:

Neuro-Immunology Laboratory University of British Columbia Hospital Room S-157, 2211 Wesbrook Mall Vancouver, B.C, Canada V6T 2B5 Centrifuge at 4 °C, transfer serum to cryovials and store at 4 °C.

Comments:

Associated with Myasthenia Gravis

This test is available exclusively to neurologists at LHSC/SJHC

Storage and Shipment:

Store serum at 4 °C and ship with ice packs. Freeze if being stored for >1 week.