

LIPOPROTEIN
(a)
SERUM/
PLASMA

Orderable - LIPOA

Turnaround Time: 10 days

Alternate Name(s):

Lpa
Lp(a)
Lipoprotein a
Lipoprotein A

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	2-10 years: 3 mL Red top Vacutainer tube
Light Green top tubes (Li-heparin) or Lavender top tubes (EDTA) are also acceptable	

Collection Information:

Minimum volume of serum/plasma required is 1 mL.

Samples must be separated within 4 hours of collection.

Interpretive Comments:

Please note change in testing method, units for results, and reference interval as of July 21, 2022.

Reference Ranges:

0.0 – 100.0 nmol/L

Storage and Shipment:

Store and send frozen.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Roche Cobas particle-
enhanced
immunoturbidimetric
assay



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)