LIPOPROTEIN (a) SERUM/ PLASMA



Pathology and Laboratory Medicine

Orderable - LIPOA

Turnaround Time: 5 days

Alternate Name(s):

Lpa Lp(a) Lipoprotein a Lipoprotein A

0

Laboratory: Core Lab

Requisition:

REQUISITION

enhanced

Test Schedule:

assay

GENERAL LABORATORY

Method of Analysis: Roche Cobas particle-

immunoturbidimetric

Referred out Monday-

Specimen:

Adult	Pediatric
5 mL Gold top	2-10 years: 3 mL Red top
Vacutainer tube	Vacutainer tube
Light Green top tubes (Li-heparin) or Lavender top tubes	
(EDTA) are also acceptable	

Collection Information:

Minimum volume of serum/plasma required is 1 mL.

Samples must be separated within 4 hours of collection.

Reference Ranges:

0.0 – 100.0 nmol/L



Thursday

Storage and Shipment:

Store and send frozen.