

PRENATAL TITRE

Orderable - PRNTITRE

Turn Around Time: 24 hours

Alternate Name(s):

Antibody Titre

Specimen:

2 x 6 mL Pink (EDTA) top Vacutainer tubes

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Reference Ranges:

See report

Interpretive Comments:

The antibody titre is a semiquantitative assay of the amount of antibody and is used primarily to monitor obstetrical patients who have produced blood group antibodies that can cause hemolytic disease of the newborn.

Comments:

Titration of clinically significant alloantibody is a reflex order by the Blood Transfusion Laboratory for pregnant women only.

Critical Information Required:

Diagnosis



Laboratory:
Blood Transfusion
Laboratory



Requisition:
[BLOOD TRANSFUSION
REQUISITION](#) or
Electronic order



Method of Analysis:
Anti-Human Globulin test



Test Schedule:
Daily
Urgent, if indicated



PRENATAL
TITRE

Storage and Shipment:

Transport to the Blood Transfusion Laboratory within 24 hours of collection at room temperature.