

BIOTINIDASE, PLASMA

Orderable - BIOTP

Specimen:

Adult	Pediatric
6 mL Green (Sodium or Lithium Heparinized) top Vacutainer tube	0-2 years: 2 x 0.5 mL Microtainers 2-10 years: 3 mL Green top tube



Laboratory:
Biochemical Genetics Lab



Requisition:
[BIOCHEMICAL GENETICS
LAB REQUISITION](#)

Collection Information:

Minimum volume of plasma required is 0.5 mL.

Store plasma at -20°C and send frozen to Biochemical Genetics.

Reference Ranges:

6-12 umol/min/L plasma



Test Schedule:
As required