



OLIGOSACCHARIDE  
SCREEN,  
URINE

## Orderable - OGCRD

---

### Specimen:

Random urine

### Collection Information:

Collect a minimum of 5 mL of urine (random) in a sterile container, then freeze before sending to lab.

Keep at -20°C

### Reference Ranges:

See report

### Storage and Shipment:

Send frozen sample to lab **immediately**.



**Laboratory:**  
Biochemical Genetics Lab



**Requisition:**  
[BIOCHEMICAL GENETICS  
LAB REQUISITION](#)



**Test Schedule:**  
Batched and tested as  
required