

OLIGOSACCHARIDE
SCREEN,
URINE

Orderable - OGCRD

Specimen:

Random urine

Collection Information:

Collect a minimum of 5 mL of urine (random) in a sterile container, then freeze before sending to lab.

Keep at -20°C

Reference Ranges:

See report

Storage and Shipment:

Send frozen sample to lab **immediately**.



Laboratory:
Biochemical Genetics Lab



Requisition:
[BIOCHEMICAL GENETICS
LAB REQUISITION](#)



Test Schedule:
Batched and tested as
required