

ORGANIC  
ACID –  
SCREENING,  
URINE

## Orderable - ORGA

---

### Specimen:

Random urine

### Collection Information:

Collect a minimum volume of 5 mL of urine (random) and send frozen to lab.

Keep at -20°C

### Reference Ranges:

See report

### Storage and Shipment:

Send frozen sample to lab **immediately**.



**Laboratory:**  
Biochemical Genetics Lab



**Requisition:**  
[BIOCHEMICAL GENETICS  
LAB REQUISITION](#)



**Test Schedule:**  
As required