

PRENATAL MICROARRAY

Orderable - Microarray aCGH

Turnaround Time: 2-3 weeks

Alternate Name(s):

Microarray
Array CGH
aCGH
Genomic Microarray

Specimen:

Amniotic Fluid
CVS
Fetal blood / Cord Blood
Maternal Blood Sample

Collection Information:

- 20 mL Amniotic Fluid
- 10-20 mg CVS (plus maternal EDTA blood 3 mL for all AF/CVS specimens)
- 2 X 1 mL fetal blood/Cord Blood (EDTA & NaHep)
- 3mL EDTA Maternal Blood Sample

Do not refrigerate, freeze or expose to heat.

Reference Ranges:

See report

Comments:

The Cytogenetics Lab is staffed from 0700-1700 (Monday-Friday), ext. 78974 (office), or 75714 (lab).

For additional information please refer to the [Cytogenetics Webpage](#).



Laboratory:
Cytogenetics Lab



Requisition:
[PRENATAL MICROARRAY REQUISITION](#) (must include patient's name, address, Ontario Health Insurance Number, PIN (if applicable), originating location, Dr's full name and address, and pertinent clinical information.

Follow up or parental samples, use the [CYTOGENETICS REQUISITION](#)



Method of Analysis:



PRENATAL MICROARRAY

Critical Information Required:

Requisition / E-order with appropriate clinical information

Storage and Shipment:

If unable to transport immediately, store at room temperature (15-25°C) and transport within 24 hours.



Test Schedule:
Batched – as required