

SKIN / FETAL TISSUE

Orderable - CGenTissue QFPCR Microarray (MD Microarray)

Turn Around Time: QF-PCR – 3 days
Microarray - 4 weeks (if clinically indicated)

Alternate Name(s):

Chromosome Analysis
QF-PCR
Microarray



Laboratory:
Cytogenetics Lab



Requisition:
[CYTOGENETICS REQUISITION](#) . (must include patient's name, address, Ontario Health Insurance Number, PIN (if applicable), originating location, Dr's full name and address, test requested, specimen type and pertinent clinical information.

[GENOMIC MICROARRAY REQUISITION](#)



Method of Analysis:

- QF-PCR
- Microarray (if clinically indicated)

Specimen:

2-3 mm² skin in a sterile container with Hank's Balanced Salt Solution (HBSS)

Collection Information:

Swab area with 70% alcohol, let dry and excise a small piece of skin (2-3 mm²) and place it in a sterile vessel with HBSS.

Reference Ranges:

See report

Special Processing:

Collect sample as aseptically as possible and send to the laboratory within 24 hours.

Critical Information Required:

Receive in laboratory within 24 hours of collection.



SKIN / FETAL TISSUE

Storage and Shipment:

If unable to transport immediately, store at 4°C and transport at room temperature (15-25°C) within 24 hours.



Test Schedule:
As required