

5-
HYDROXINDOLE
ACETIC
ACID,
URINE

Orderable - 5HIAAU24

Turnaround Time: 5 days

Alternate Name(s):

5-HIAA

Specimen:

24-hour urine

Collection Information:

[24-Hour Urine Diet and Medication Restrictions: 5-HIAA, Catecholamines, Metanephrines, VMA/HVA, Serotonin](#)

Collection of 24-Hour Urine Sample

Collect urine over a 24-hour period in a container with 30 mL of 6 N HCl for adult samples or 15 mL of 6 N HCl for pediatric samples. The final pH must be maintained at 2-4.

Reference Ranges:

≤ 43.0 μmol/day

Special Processing:

Transfer an aliquot of the measured urine to a 5-mL Greiner tube.

Storage and Shipment:

Store and ship urine aliquot refrigerated.



Laboratory:
Toxicology/Special
Chemistry Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
LC-MS/MS



Test Schedule:
Monday – Friday
0800-1600