

# SIROLIMUS, WHOLE BLOOD

## Orderable - SIR

Turnaround Time: 1 day

STAT: Same day if received in Core Lab by 1000 h

### Alternate Name(s):

Rapamycin

### Specimen:

Adult	Pediatric
4 mL K <sub>2</sub> or K <sub>3</sub> EDTA Lavender top Vacutainer tube	0-2 years: 0.5 mL Lavender Microtainer 2-10 years: 3 mL Lavender top

### Collection Information:

Send to Lab immediately after collection.

Clotted blood not suitable for analysis.

### Reference Ranges:

No established reference range. Concentrations are measured in ug/L/

### Storage and Shipment:

Ship between campuses on ice packs.

Referred in samples: Ship on ice packs.



**Laboratory:**  
Toxicology/Special  
Chemistry Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
LC-MS/MS



**Test Schedule:**  
Monday--Saturday  
Specimens received in  
Core lab after 1000 h will  
be processed the next  
working day